

# PRODUCT INFORMATION

## FFP3 NR D preformed respiratory mask with valve by Uvex



### Scope of application and properties

- + **Maximum protection and wearing comfort:** Particle-filtering cup shaped mask with exhalation valve. The FFP3 respiratory mask provides a high-level protection against particles, aerosols and pollutants, as well as protection against airborne infectious agents with a retention capacity of more than 99%. The exhalation valve ensures minimal breathing resistance (even in the presence of heavy dust) and reduces heat and moisture build-up under the mask. The integrated nose clip can be flexibly adjusted, the 4-point headband attachment, the soft material edges and the all-round sealing lip ensure a reliable tight fit and wearing comfort. The FFP3 respiratory mask is available in 2 sizes. Type tested and certified in accordance with PPE Regulation (EU) 2016/425 (Personal protective equipment).
- + **Application area:** Protection from solid and liquid particles, harmful or non-toxic particles, e.g. CMR-drugs<sup>1)</sup> (also protection from biological agents).
- + **Protection capacity:** Checking the fit of your mask can be easily by carrying out a simple test (the filter material moves when breathing in). No guarantee for all hazardous substances, CMR-drugs or chemicals. Does not protect against gases and vapours.
- + **Replacement interval:** Depends on inhalation resistance and personal well-being. Maximum use for one shift (8 hours) with at least three breaks in between. For use with CMR-drugs, a maximum of 4 hours, as long as there is no obvious contamination. Replace immediately in the event of visible contamination. Single use – Do not re-use respiratory mask!
- + **Before use:** Check for any damage. Do not use damaged masks.
- + **Limitation:** Single use mask! Laboratory test results can differ from values obtained in practice considerably and, as a result, the useful life of the mask may increase or decrease. Facial hair and scars can impair the fit of the mask and therefore the level of protection.
- + **Disposal:** If used in combination with CMR-drugs and contaminated areas/items: waste requiring supervision (waste code: 18 01 04 in accordance with 2000/532/EC); in the event of heavy contamination, waste requiring special supervision (waste code: 18 01 08\* in accordance with 2000/532/EC); collect and dispose of waste separately!  
1) Carcinogenic; mutagenic; toxic to reproduction

### Type

Item No.	121275	121276
Size	Universal (adapted for larger faces)	S (adapted for smaller faces)
Packaging unit (PU)	1 mask	

### Material Properties

Latex free	Yes, complete mask
Filter/mask	Melt blown nonwoven polypropylene, polyester
Sealing lip	PVC (Polyvinyl chloride) foam
Valve	ABS (Acrylonitrile butadiene styrene)
Headband	Textile

---

**Nose clip**

Aluminium

---

**Weight**

Excl. Packaging: 29 g

---

## Filter Properties

---

Retaining capacity

At least 99%

---

Dolomite dust test requirements FFP3  
respiratory masks

Passed

---

Protection against dust, smoke and  
aerosols in water or oil

Up to 30-fold of threshold

---

Gases

Unsuitable

---

## CE-Marking

CE-mark for complex PPE in category III in accordance with PPE Regulation (EU) 2016/425. The type test was based on EN 149:2001+A1:2009. Documented by the EC type test certificate of the manufacturer: Control measures by intermediary notified body BSI Group The Netherlands B.V. (2797).

---

## Notified Body „2797“

BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands.

---

## Quality Management System

Our quality management system is tested and certified by TÜV Product Service GmbH (a certification body accredited by the German Accreditation Council) in accordance with DIN EN ISO 9001:2015. Regular audits and production inspections guarantee the quality of our products.

---

## Storage and Transport Conditions

- + Dark (protection from direct UV light and sunlight)
  - + Cool, dry (Relative humidity 30% - 60%)
  - + No contact with pointed and/or sharp objects
- 

## Shelf Life

3 years

---

## Manufacturer

Uvex Arbeitsschutz GmbH, Würzburger Str. 181 – 189, 90766 Fürth, [www.uvex-safety.com](http://www.uvex-safety.com)

---

## Distributor

Berner International GmbH, Werner-von-Siemens-Str. 19, 25337 Elmshorn, Tel: +49 4121 43560, Fax: +49 4121 435620  
[info@berner-safety.de](mailto:info@berner-safety.de), [www.berner-safety.de](http://www.berner-safety.de)

---