

PRODUCT INFORMATION

ISYSOFT Cleaning instructions

Cleaning guide for CMR drugs contaminations (e.g. cytostatics)

For the daily cleaning of safety cabinets the BERNER ISYSOFT cleaning cloths should be used in 2 steps:

Step 1: Clean the surface with an alkaline – cleaning solution > pH 11 (e.g.: sodium hydroxide solution (0.1 mol), Mucosol (2%-solution), etc.

Step 2: Subsequently clean with Isopropanol (IPA) 70%.

Results of the investigations of IUTA e.V. as part of the validation measures: Degree of cleanliness > 99.95%

Cleaning of highly contaminated surfaces surfaces with cytostatics (e.g. after accidents)*:

For highly contaminated surfaces please repeat 3 times the Steps 1 and 2:

Step 1: Clean with an alkaline cleaning solution > pH 11 (e.g.: sodium hydroxide solution, Mucosol,)

Step 2: Clean with Isopropanol (IPA) 70%

Step 1: Clean with an alkaline cleaning solution > pH 11 (e.g.: sodium hydroxide solution, Mucosol,)

Step 2: Clean with Isopropanol (IPA) 70%

Step 1: Clean with an alkaline cleaning solution > pH 11 (e.g.: sodium hydroxide solution, Mucosol,)

Step 2: Clean with Isopropanol (IPA) 70%

Results according to the investigation of the IUTA e. V. within the scope of validation measures how:

Cleaning results in conjunction with the BERNER Isysoft cleaning pads: > 99.95%.

Caution: This manual is not suitable for flooring made of natural materials (eg linoleum). In this case, use an alkaline cleaner <pH 10. * = This cleaning method is also recommended by the DGOP (see current QUAPOS).