

# PRODUCT INFORMATION

## FFP-3 Respirator Mask Dräger X-plore 1730



Dräger X-plore 1730 FFP3 NR D (without exhalation valve)

### Scope & properties

- + **Features:** Filtering folding half-mask without exhalation valve. The soft material and special shape with flexible nose clip ensures that it can be easily adjusted to fit any face shape. Comfortable to wear. The threaded elastic strap makes it easy to put on and take off. Type tested and certified in accordance with the EU PPE regulation (EU) 2016/425 (personal protective equipment) for PPE category III.
- + **Area of application:** Protection from solid and liquid particles, e.g. for handling CMR<sup>1)</sup> drugs (e.g. cytostatics, virostatics, biological agents and viruses).
- + **Protection capacity:** Protection from all CMR drugs or chemicals is not guaranteed. Does not protect against gases and vapours!
- + **Replacement interval:** Depends on type of activity and individual comfort. Maximum use for one shift (8 hours) with at least three breaks in between. For use with CMR-drugs, a maximum of 4 hrs, so far as there is no obvious contamination. Replace immediately in the event of visible contamination. Single use mask – Do not re-use respirator mask!
- + **Before use:** Check for damage. Do not use damaged masks.
- + **Limitation:** Single use mask! Laboratory test results can differ from values obtained in practice considerably and, as a result, the useful life of the mask may increase or decrease. Facial hair and scars can impair the fit of the mask and therefore the level of protection.
- + **Disposal:** If used in combination with CMR-drugs and contaminated areas/items: waste requiring supervision (waste code: 18 01 04 in accordance with 2000/532/EC); in the event of heavy contamination, waste requiring special supervision (waste code: 18 01 08\* in accordance with 2000/532/EC); collect and dispose of waste separately!

<sup>1)</sup> Carcinogenic; mutagenic; toxic to reproduction.

### Type

Size	Universal
Item No. (non-sterile)	121289
Packaging unit (PU)	1 piece
Size of single pouch (PU)	155 mm x 115 mm

### Material Properties

Design	Several filter material layers, non-woven
Nose clip	Tin plate, aluminium-free
Straps	Latex straps, with textile covering
Weight	Excl. packaging : 14.3 g

## Filtration efficiency

Test aerosol	Minimum filtration efficiency
Sodium chloride	> 99%
Paraffin oil (as per EN 149)	> 99%
Dolomite dust test	Passed
Gases	Unsuitable

## CE-Mark

CE-mark for complex PPE in category III in accordance with PPE Regulation (EU) 2016/425. The type test performed was based on EN 149:2001+A1:2009. Documented by the EC type test certificate of the manufacturer. Quality assurance (EC quality assurance system with monitoring) by notified body DEKRA EXAM GmbH (0158) in accordance with (EU) 2016/425.

## Notified Body 0158

DEKRA EXAM GmbH, Dinnendahlstraße 9, 44809 Bochum, Germany

## Quality Management System

Our quality management system is tested and certified by TÜV Product Service GmbH (a certification body accredited by the German Accreditation Council) in accordance with DIN EN ISO 9001:2015. Regular audits and production site inspections guarantee the quality of our products.

## Storage and Transport Conditions

- Dark (protect from direct UV light and sunlight)
- Cool
- Dry (Relative humidity 30% - 60%)
- No contact with pointed and/or sharp objects

## Shelf Life

4 years

## Manufacturer

Dräger Safety AG & Co. KGaA, Revalstr. 1, D-23560 Lübeck

## Distributor

**Berner International GmbH**, Werner-von-Siemens-Str. 19, D-25337 Elmshorn, Tel: +49 4121 43560  
[info@berner-safety.de](mailto:info@berner-safety.de), [www.berner-safety.de](http://www.berner-safety.de)