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Referenz

1. International Council for Harmonisation (ICH). ICH Q9 - Quality risk management. In:2005.
2. International Council for Harmonisation (ICH). PHARMACEUTICAL DEVELOPMENT Q8(R2). In: INTERNATIONAL CONFERENCE ON HARMONISATION OF TECHNICAL REQUIREMENTS FOR REGISTRATION OF PHARMACEUTICALS FOR HUMAN USE, ed2009.
3. Pavičić M, Wagner T. Risk & Science-Based Validation of Cleanroom Garments (<http://www.ivtnetwork.com/article/risk-science-based-validation-cleanroom-garments>). Journal of GXP Compliance (IVT Network). 2019;23(4).
4. European Medicine Agency (EMA). Annex 1 - Manufacture of Sterile Medicinal Products - draft for public consultation. In: European Medicine Agency (EMA), ed. Brussels: Directorate-General for Health and Food Safety, Unit SANTE B/4, BE-1049 Brussels; 2017.

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T. +352 3666 5111
mycustomerservice.emea@dupont.com

www.ipp.dupont.com

Risk and science-based validation of cleanroom garments

Milenko Pavicic and Thierry Wagner

Abstract

Important quality attributes of cleanroom garments that to worn during the manufacture of sterile medicines include: cleanliness; sterility; particle and microbial filtration efficiency; durability; usability; and comfort. Important risk factors related to cleanroom garment systems include gowning processes and related activities, such as laundering, packing, sterilization, repairs, storage, handling and logistics, as well as change management. Because many factors contribute to the overall quality, adequacy and suitability of cleanroom garment systems, a thorough and focused validation of cleanroom garments is critically important.

After providing a review of current and emerging regulations and standards, this article proposes a risk- and science-based quality-by-design approach for the development, implementation and validation of sterile cleanroom garment systems. With this approach, more effort is spent at the front-end during the design phase as well as during design qualification. This will lead to designed-in risk reductions; enhanced scientific knowledge on selected technical solutions; and better awareness of limitations and residual risks. As a result, there should be fewer issues during cleanroom qualifications and process validations, leading to more effective routine operations as well as improved patient safety. The proposed approach, if implemented correctly, is not only the correct strategy to effectively control contamination risks related to people, but also an adequate response to the latest regulatory requirements.

Introduction

Parenteral medicines are administered through injection, infusion or implantation, and must be sterile and pure to assure product safety. Therefore, the manufacture of parenteral medicines requires a controlled and validated clean production environment. Sterile medicines can be manufactured by terminal sterilization. In this case,

the product is sterilized in its final packaging, resulting in a sterility assurance level (SAL) of at least 10^{-6} . If terminal sterilization of the final product is not possible, aseptic manufacturing is the only alternative. In aseptic production, exposure of the sterile product to the environment may take place during different stages of the manufacturing process (i.e., from preparation of the product using sterile ingredients to filling the product in its final container). In Europe, exposure of sterile products to the environment is only allowed in EU GMP grade A zones placed in a grade B cleanroom or in an isolator.¹

Aseptic manufacturing of sterile products requires a high level of contamination control. Figure 1 shows the important elements that must be addressed in a contamination control strategy to assure purity of the

products, sterility of sterile products and good microbiological quality of non-sterile products.

According to current Good Manufacturing Practices (GMP) guidelines, (EU,² US,³ Japan,⁴ WHO^{5,6}), processes, equipment, facilities and manufacturing activities should be managed in accordance with Quality Risk Management (QRM) principles⁷ that provide a proactive means of identifying, scientifically evaluating and controlling potential risks to quality. A quality-by-design approach,⁸ combined with effective risk management, should be applied to ensure the safety, quality and efficacy of sterile products. This comes from the belief that quality cannot be tested into the product, it can only be built into the design of products and the processes used to produce them. Risk assessments must be performed to identify, assess, eliminate and control


Figure 1: Important elements of a contamination control strategy for sterile manufacturing

contamination risks from the design phase of an aseptic process; to monitor and detect contamination; and to establish process requirements and acceptance criteria for all elements of a sterile manufacturing process. Risk assessments should be documented as well as maintained and should include the rationale for decisions taken in relation to mitigating risks, discounting of potential risks and residual risk.

An important risk factor in sterile manufacturing is personnel. People can cause contamination of the production environment and products in many ways.⁹ Contamination from people mainly consists of hair, skin cells, skin flakes, saliva, sebaceous matter, sweat, particles from clothing and many different exogenous particles and substances picked up in the environment. These contamination sources mostly contain different endogenous (i.e., commensal) and exogenous microorganisms present in numbers which vary from a few (e.g., skin cell) to thousands (e.g., sweat) or even millions (e.g., saliva). Therefore, it is important that people working in an environment, where sterile products are manufactured, wear adequate cleanroom garments.

Contamination control related to people starts with good personal hygiene based on adequate hygienic procedures, such as hand washing and hand disinfection; aseptic behavior; aseptic skills; and a good working discipline. Adequate cleanroom garments, as well as undergarments are critically important to minimize the risk of contaminating the environment or products with contamination generated by people. Cleanroom undergarments¹⁰ serve as a first barrier against contamination from people. The cleanroom garments must form a very robust barrier between the person and the environment.

Additional protective measures, such as adequate cleanroom footwear, cleanroom socks, head coverings, face masks, eye coverings and (sterile) gloves, may be necessary or required from a regulatory point of view, to minimize the contamination risk.

Important quality attributes of cleanroom garments are cleanliness (free from chemicals, particles, pyrogens, fibers); sterility; particulate and microbial filtration efficiency; durability (tear, puncture, seam strength, abrasion); and comfort.

Personal protection against chemical or biological agents may also be a relevant quality attribute. A quality-by-design approach⁸ combined with effective risk management,⁷ should also be applied to the design, selection and implementation of adequate cleanroom garments.

With the quality-by-design approach, more effort is spent at the front-end, in the design phase as well as in the design qualification, leading to designed-in risk reductions; better understanding of key aspects, limitations and residual risks; and fewer issues during final simulation runs and routine operations. This approach also creates the basis for proper root cause analysis (in case of issues) and adequate change management.

Other important risk factors related to cleanroom garments include gowning procedures and processes and activities related to cleanroom garments, such as laundering, packing, sterilization, repairs, storage, handling and logistics. Because many factors contribute to the overall quality, adequacy and suitability of cleanroom garment systems, a risk- and science-based validation of cleanroom garment systems is very important.

This article provides an overview of the various qualification, validation and monitoring aspects of cleanroom garments.

Current regulatory guidance for validation of cleanroom garments

Depending on the jurisdiction, aseptic production of sterile medicines must meet various regulatory requirements such as those set out in Annex 1 of the EU Guidelines to Good Manufacturing Practice¹ or in the U. S. Food and Drug Administration (FDA) Guidance for Industry on sterile drug production³.

The FDA guidance outlines various recommendations for gowns, such as proper gown control; no unreasonable contamination risk to the gown; providing a barrier between the body and exposed sterilized materials; and preventing contamination from particles generated by, and microorganisms shed from the body. Gowns should be sterilized and non-shedding. Methods used to don each gown component in an aseptic manner should be detailed. Manufacturers should implement an aseptic gowning qualification program to assess the ability of a cleanroom operator to maintain the quality of the gown after performance of gowning

procedures, including periodic requalification and microbiological monitoring of strategically selected locations on the gown.

The current EU GMP guidelines require sterile (sterilized or adequately sanitized) garments to be provided for grade A/B areas and changing and washing to follow a written procedure designed to minimize contamination of clean area clothing or carry-through of contaminants to the clean areas. It further requires that clothing and its quality be “appropriate” without defining what would be considered appropriate. It also requires clothing to “be worn in such a way as to protect the product from contamination”. In terms of attributes, “protective clothing should shed virtually no fibers or particulate matter and retain particles shed by the body”. Reusable garments are required to be cleaned and handled in such a way that the garment does not gather additional contaminants that can be shed later.

The current EU GMP Annex 1 for the Manufacture of Sterile Medicinal Products¹¹ includes little guidance on cleanroom garment qualification except that it needs to be “appropriate”.

The new draft EU GMP Annex 1¹² published for consultation in December 2017 explicitly introduces the application of QRM principles and provides more details on gowning, including the requirement that gowning is part of a holistic contamination control strategy. It further requires that:

- Personnel are trained on gowning practices, which must be assessed
- Personnel are qualified through a successful aseptic process simulation test
- Microbial monitoring of personnel is performed

In addition, this new draft EU GMP Annex 1¹² requires that garments are visually checked for cleanliness and integrity. It also provides several new requirements regarding clothing of grade A/B; body parts that should be covered; the prevention of droplets, particles and fibers; and the folding of packaged garments to make sure that “contact to the outer surfaces of the garment is reduced to a minimum”.

A key addition is the requirement that “reusable garments should be

Main feature

replaced at a set frequency determined by qualification or if damage is identified". This requires manufacturers to produce data regarding the effect of reprocessing on the fabric and the overall garments.

However, there is little guidance on how to qualify those garments other than stating in subclause 4.11 that "clothing and its quality should be appropriate for the process and the grade of the working area. It should be worn in such a way as to protect the product from contamination".

The Japanese Guidance on the Manufacture of Sterile Pharmaceutical Products by Aseptic Processing,⁴ developed by a task force of Japanese experts, is an extensive document covering all areas of aseptic processing. It includes recommendations on gowning, operations, monitoring and controls that are comparable to other standards. It provides some interesting design recommendations on gowning and de-gowning areas and highlights the need to establish appropriate control procedures, including visual inspection. It also defines maximum allowable frequency of steam sterilization for reusable materials, such as aseptic gowns, to "ensure maintenance of

specifications, safety, and intended functions after repeated exposure to steam at its maximum intensity".

ISO 14644-5¹³ includes an informative annex B on cleanroom clothing requirements that provides some useful guidance on aspects to consider during qualification of such clothing. It provides guidance on barrier properties; evaluation of electrostatic properties; some guidance on fit and function; and helpful guidance on construction, finishing of seams and general design criteria that can be used to establish URS. It proposes use of the dispersal chamber or body box¹⁴ as a simulation procedure to evaluate performance; recommends that shelf life of sterile packaging be determined; and includes considerations on thermal comfort and some guidance on cleaning, referring to IEST-RP-CC003.4.¹⁴

ISO 13408-1¹⁵ on general requirements of aseptic processing includes some general requirements on cleanroom garments but does not provide much guidance on cleanroom garment system qualifications.

IEST-RP-CC003.4 "Garment Systems Considerations for Cleanrooms and Other Controlled Environments"¹⁴ provides guidance on design, selection,

specification, maintenance and testing of garment systems. It includes useful guidance on material qualification attributes, considerations on processing (cleaning, re-sterilization, etc.), gowning system specifications and quality management. Appendix B proposes various tests for assessments of particle penetration and garment cleanliness, which includes the well-known and very useful body box test,^{16,17} that allows simulation of the actual use of the garment, as well as the Helmke drum test.¹⁸ Overall, IEST-RP-CC003.4¹⁴ is the most useful document to support qualifications of cleanroom garment systems.

The EU general guidance on validation (GMP Annex 15¹⁹) provides the general framework that we will apply to the qualification of cleanroom garment systems.

Validation approach for cleanroom garments

The main stages of validation of equipment, facilities, utilities or systems are:

- Definition of User Requirements Specification (URS)
- Design Qualification (DQ)


Figure 2: Overview of the different validation stages for cleanroom garments used in EU GMP grade A/B cleanrooms

- Installation Qualification (IQ)
- Operational Qualification (OQ)
- Performance Qualification (PQ)

This approach is also appropriate for cleanroom garments. Some stages of the validation can focus mainly on the quality of the cleanroom garment itself, but in other stages, the other items of the cleanroom clothing system (i.e., cleanroom undergarments, footwear, socks, head coverings, face masks, eye coverings, gloves and other accessories) must be included. The packaging (sterile and non-sterile barriers) of the cleanroom garments should be part of the validation. An overview of the different validation stages and validation items for cleanroom garments is given in Figure 2. Each validation stage must be formally finalized before progressing to the next stage.

The GMP (EU, US, Japan) guidelines state that QRM7 should be used for qualification and validation activities. If changes occur during the project phase or during commercial production (e.g., change of garment design, fabric, zipper type, packaging, laundering process or sterilization process), the risk assessments must be repeated to determine if additional validation must be performed. The way in which risk assessments are used to support qualification and validation activities should be clearly documented.

For critical goods and materials, such as cleanroom garments, it is also important to qualify the supplier. This qualification should provide an appropriate level of confidence that the supplier is able to supply cleanroom garments with consistent quality and acts in compliance with regulatory requirements. A supplier qualification should also include qualification of subcontractors, suppliers of base materials (e.g., fabric and garment accessories) and outsourced service providers (e.g., laundries and sterilization facilities).

User Requirements Specification (URS)

The specification for cleanroom garments should be defined in a User Requirements Specification (URS). The URS is a document that specifies requirements necessary and sufficient to create a feasible design, meeting the intended purpose of the material,

equipment, utility or system. The URS may also include additional requirements, such as protection of people against chemical and/or biological agents.

An example of a URS for cleanroom garments is given in Table 1.

Cleanroom garments for use in EU GMP grade A/B cleanrooms should

Table 1: Example of a URS for cleanroom garments used in EU GMP grade A/B cleanrooms

| General |
|---|
| 1. The clothing must consist of a suit or coverall; a hood that can be firmly and reliably tucked into the neck of the suit; and boots. The hood may also be attached to the suit. It must be possible to tuck the trouser legs firmly and reliably into the boots. Boots shall have an anti-slip sole. |
| 2. The garments may be reusable or single-use. |
| 3. Unpacking and gowning must be conducted in an aseptic manner. |
| 4. The garments must fit well and be available in different sizes (e.g. XS to XXL). |
| 5. The clothing must be comfortable to wear (e.g., ease of movement, thermal, tactility). |
| 6. Garments must be traceable, including number of laundering and sterilization cycles. |
| Cleanliness and sterility |
| 1. The garments should be clean (low levels of ions and organic extractables). |
| 2. The garments shall not shed fibers, including filaments of the fabric and thread. |
| 3. The garments shall not shed too many particles ($\geq 0.5\mu\text{m}$ and $\geq 5\mu\text{m}$). Maximum number of particles ($\geq 0.5\mu\text{m}$ and $\geq 5\mu\text{m}$) must be specified. |
| 4. The garments shall be sterile (sterility assurance level $\leq 10^{-6}$). |
| 5. The garments shall have a low pyrogenicity. |
| Barrier properties |
| 1. The garments shall have a good particulate and bacterial filtration efficiency. |
| Packaging |
| 1. The garments shall be double or triple sterile packed. |
| 2. The sterile packaging must be robust and should endure common manipulations without jeopardizing the integrity of the packaging. |
| 3. The sterile packaging shall be clean. |
| 4. The sterile packaging shall not shed too many particles ($\geq 0.5\mu\text{m}$ and $\geq 5\mu\text{m}$) and shall not shed fibers. Maximum number of particles ($\geq 0.5\mu\text{m}$ and $\geq 5\mu\text{m}$) must be specified. |
| 5. The space between the sterile packaging shall not contain too many particles ($\geq 0.5\mu\text{m}$ and $\geq 5\mu\text{m}$) and should not contain fibers. Maximum number of particles ($\geq 0.5\mu\text{m}$ and $\geq 5\mu\text{m}$) must be specified. |
| 6. The sterile packed garments must be packed into a clean, low bioburden bag that is present in a firmly closed clean box. |
| 7. The box and the last sterile packaging of the garments must be labeled with a clear description of the manufacturer, content, size, batch number, production date, expiry date, storage conditions and indicator that the sterilization was done. |
| Identification |
| 1. Each garment must have a unique identification number. |
| 2. Labels must be firmly attached to the packaging or exist as an integral part of the packaging. |
| 3. Labels must be clean and should not shed too many particles ($\geq 0.5\mu\text{m}$ and $\geq 5\mu\text{m}$) and should shed no fibers. Maximum number of particles ($\geq 0.5\mu\text{m}$ and $\geq 5\mu\text{m}$) must be specified. |
| 4. Labels and printed text must be resistant to disinfectants that are commonly used (e.g., 70% alcohol, 6% hydrogen peroxide solution or a 600-ppm hydrogen peroxide vapor, 1000 ppm active chlorine). |

adequately protect the environment and products against contamination from people. A trained and qualified operator, wearing a nonwoven polypropylene head cover, clean polyester two-piece undergarment, clean dedicated socks, double sterile gloves, a sterile face mask and sterile goggles, must be able to work at least three hours in the same set of cleanroom garments without causing unacceptable (cGMP) levels of contamination of the garments and the aseptic working environment.

Design Qualification (DQ)

During DQ, compliance of the cleanroom garment design with cGMP must be demonstrated and documented, and the requirements of the URS must be verified. The purpose of the DQ is to confirm that the selected cleanroom garment is qualified for the intended use. Therefore, it should also include tests to simulate the intended use. The DQ must be executed and authorized by suitably qualified persons who are knowledgeable enough to challenge the proposed design and its performance.

Following the model of design validations of sterile barrier systems described in ISO 11607-1,²⁰ it is recommended to split the DQ into four key areas: material qualification, performance testing, stability testing and usability evaluation. For reusable

garments, reprocessing should be the subject of a separate DQ by the manufacturer and IQ-OQ-PQ by the supplier.

The material qualification includes the qualification of key characteristics and key properties of the materials and fabrics used, the cleanroom garments and the packaging.

Performance testing includes testing of the cleanroom garments and the packaging under simulated and standardized conditions using standardized test methods.

Stability testing must be performed to assure that key material characteristics and properties remain sufficiently stable during the life cycle. Characteristics and properties that change over time (e.g., deterioration of filtration efficiency of garments due to repeated laundering and sterilization; wear of the garments due to multiple use; and changes to the integrity of the sterile packaging during storage due to long-term effects of gamma irradiation) should be validated under worst-case conditions. Information for material qualification, performance testing and stability testing normally should be provided by the supplier. It is important to verify that the data has been generated using validated and sound scientific methods.

The purpose of the usability evaluation is to assure that the cleanroom

garments can be used with acceptable remaining contamination and safety risks due to the design of the garments and the established gowning, working and de-gowning procedures. The usability evaluation is typically done by the end-user. However suppliers can also evaluate their garments for the intended use and supply that data to users for verification and further mitigation of identified risks during gowning and operations. The concept of usability engineering and testing has developed into a well-accepted way to successfully reduce use-related risks by reviewing these risks during the design phase and systematically reviewing the design and use of the product against those identified risks (see also IEC 62366-1 on usability engineering of medical devices²¹).

Relevant items to be covered for each of the four areas are presented in Table 2. A summary is given in Figure 2.

Reusable versus single-use cleanroom garments.

The validation of reusable cleanroom garments is more complex and more extensive compared to single-use cleanroom garments. Repeated laundering, repeated sterilization, multiple use and repairs influence the quality of reusable cleanroom garments. This also means that the influence of these factors must be validated throughout the entire

Table 2: The four key areas of the Design Qualification (DQ) for cleanroom garments used in EU GMP grade A/B cleanrooms

| Material qualification | Performance testing | Stability testing | Usability evaluation |
|---|--|--|--|
| Cleanroom garments <ul style="list-style-type: none"> • Fiber and particle shedding • Sterilization compatibility • Sterility assurance level • Pyrogenicity • Particle filtration efficiency • Bacterial filtration efficiency • Porosity • Surface resistivity • Perforation resistance • Mechanical strength • Chemical resistance • Protection against biological agents | Cleanroom garments <ul style="list-style-type: none"> • Body box testing • Helmke drum test | Single-Use garments <ul style="list-style-type: none"> • Properties and characteristics at the end of shelf life Reusable garments <ul style="list-style-type: none"> • Properties and characteristics after maximum number of laundering and sterilization cycles | Use scenarios <ul style="list-style-type: none"> • Transfer to classified storage area • Readability of label • Easy opening of packaging • Aseptic unfolding of garments • Gowning • Donning additional accessories (e.g., sterile gloves, face mask, goggles) • Work situations • Safety, biosafety • De-gowning |
| Packaging <ul style="list-style-type: none"> • Fiber and particle shedding • Bioburden • Penetration of commonly used disinfectants | Sterile packaging <ul style="list-style-type: none"> • Influence of transport on integrity/sterility (ISO 11607-1) | Sterile packaging <ul style="list-style-type: none"> • Packaging integrity/sterility at the end of shelf life (ISO 11607-1) | Packaging <ul style="list-style-type: none"> • Aseptic presentation of garments (multiple layers) |
| Sterile packaging <ul style="list-style-type: none"> • ISO 11607-1 | | | |

life cycle in the frame of the stability testing and the performance testing at the end of the life cycle. In addition, not only the supplier of the garments but also the cleanroom laundry, sterilization facilities and repair service must be qualified. Reprocessing should be the subject of a separate DQ by the manufacturer and IQ-OQ-PQ by the supplier.

Installation Qualification (IQ)

The IQ for cleanroom garments is a formal check to verify if all required elements of the cleanroom gowning system are present. These include the gowning and de-gowning facilities; certificates of conformance and/or analysis; implementation of instructions from the supplier; standard operating procedures for gowning and de-gowning; logistical processes for the cleanroom garments and related accessories; and the operator training and qualification plan. Risk assessments that were executed as part of the DQ of the cleanroom garments should be finalized and risk controls should be implemented.

In addition, it is important to verify that the correct materials have been received for performing the OQ and PQ (i.e., the correct cleanroom garments, correctly folded, in the correct packaging and correctly labeled). A summary of items to be included in the IQ is given in Figure 2.

Operational Qualification (OQ)

During the OQ, the objective is to qualify the gowning and de-gowning concept. For this purpose, all relevant steps of the gowning and de-gowning process, including logistics, should be qualified. In addition, the aseptic presentation of the garments should be qualified. To validate the aseptic gowning procedure, at least three independent, consecutive, successful visual and microbiological assessments for at least one person who is trained for aseptic gowning should be performed. The OQ should also include a formal assessment to verify that different work tasks can be executed properly from a practical point of view (e.g., moving, bending, stretching and lifting). It is recommended to perform this assessment with all available sizes of the cleanroom garments and with people of different body shapes. A summary of items to be included in the OQ is given in Figure 2.

Performance Qualification (PQ)

During PQ, the objective is to validate the performance of the cleanroom garment system when it is actually in use. The requirements specified in the URS must be complied with fully.

The PQ of cleanroom garments consists of two stages. In the first stage of the PQ, compliance with aseptic gowning procedures should be assessed and confirmed. This gowning qualification must involve both a visual and microbiological assessment. The visual assessment is to qualify that people don the cleanroom garments in a correct and aseptic manner, which shall be described in detail in a gowning procedure. After gowning, the microbiological quality shall be assessed by taking surface samples from several locations on the cleanroom garments, gloves, goggles and face mask. Locations must be determined based on a risk assessment.

Each person accessing an EU GMP grade A/B environment must perform a gowning qualification. For the PQ, it is important to determine how many gowning qualifications are required to demonstrate compliance with the requirements. Typically, initial gowning qualification is performed three times for each person. It is important that adequately trained, qualified and experienced persons execute the PQ to exclude failures due to causes other than quality issues with the cleanroom garments.

The second stage of the PQ focuses on the validation of the microbiological quality of the gowned personnel with the garments and other accessories (e.g., gloves, face mask, goggles) during the actual work (e.g., aseptic compounding, aseptic filling, cleaning and disinfection, and other activities).

The second stage of the PQ also includes validation of the microbiological and particulate quality of the environment people are working in and the execution of aseptic process validations (i.e., media simulations or media fills). The number of runs for these validations must be determined based on a risk assessment. Typically, these validations are performed three times. To exclude failures due to causes other than quality issues with the cleanroom garments, it is important that adequately trained, qualified and experienced personnel execute the PQ in areas with an excellent quality history.

The PQ is typically done under worst-case conditions. These worst-case conditions must be determined based on a risk assessment. Also, the actions that should be taken if established criteria are not met during the PQ, must be defined before executing the PQ. Only after a successful PQ can the cleanroom clothing system be formally implemented. A summary of items to be included in the PQ is given in Figure 2.

Revalidation and change management

The cleanroom clothing system should be evaluated at an appropriate frequency (e.g., annually or biennially) to confirm that it remains in a state of control. Gowning qualifications shall be repeated at least annually and even more frequently in cases where there is doubt about the quality of the aseptic gowning process or aseptic gowning skills of specific persons. The cleanroom clothing system is included in validations which must be performed periodically (i.e., cleanroom qualifications under dynamic conditions and aseptic process validations).

Changes must be reviewed critically and may lead to revalidations that are more or less extensive, depending on the type of change. Properly and well documented DQs, as well as IQ-OQ-PQ, are the basis for successful change management.

Monitoring

Personnel monitoring must be part of the environmental monitoring program^{1,3,22}. The microbiological quality of cleanroom garments for persons working in a grade A/B environment must comply with the EU GMP1 grade B limit for surface samples (i.e., the action limit is 5 CFU/contact plate). Alert limits are usually lower (e.g., 2 or 3 CFU/contact plate). Cleanroom garments, face masks, goggles and gloves are typically sampled at the conclusion of activities in a grade A/B area, but just before leaving the area. For this "exit monitoring," contact samples are taken from different locations. Sample locations must be determined based on a risk assessment. After sampling, the person must leave the area to prevent spreading contamination due to medium residues present on the cleanroom garments.

Gloves should be sampled after performing activities in a grade A environment to verify the quality of the aseptic conditions and aseptic handling. Gloves of operators working in a grade B environment should also be monitored during each work shift. Gloves are typically sampled with a frequency ranging from once to multiple times per work shift.

In addition to personnel monitoring, samples from several locations in the cleanroom should be taken to determine if the production environment and processes are in control. Sampling is typically performed during (passive and active microbiological air sampling and particle counting) or at the conclusion (surface sampling) of operations but may also be performed under static conditions (i.e., in the at-rest state of an area), to verify cleanliness.

In the case of non-conformities, assessments must be done to determine root causes. The cleanroom clothing system as a potential root cause should be included in these assessments.

It is also recommended to assess if gowning procedures, cleanroom behavior guidelines and aseptic procedures are followed correctly. These visual assessments should be done on a regular basis.

Conclusion

A science- and risk-based quality-by-design approach for the development, implementation and validation of sterile garment systems for EU GMP grade A/B aseptic processing areas is not only the correct approach to effectively control contamination risks related to people, but also an adequate response to the latest regulatory requirements. The new EU GMP Annex 1 draft is based on QRM principles and introduces the concept of a holistic contamination control strategy that considers all aspects of contamination control over the entire life cycle based on thorough technical knowledge and sound process know-how. Considerable efforts will be required by manufacturers to update their technical files, with cleanroom garments being just one of the many aspects.

With a risk-based quality-by-design approach applied to cleanroom garment systems, more effort is spent at the front-end, in the design phase, as well as in the design qualification. This will lead to designed-in risk reductions;

better scientific knowledge of key aspects, attributes, limitations and residual risks of the selected technical solutions; and fewer issues during cleanroom qualifications, process validations and routine operations. In case of failures, it can be difficult to determine the root cause or the elements that have failed. That's why a quality-by-design approach, with focused and extensive design qualifications for each element, is the only way to successfully and systematically reduce the risk of failure. This approach also creates the foundation for adequate and risk-based change management.

References

1. European Commission. Eudralex - Volume 4 - EU Guidelines to Good Manufacturing Practice - Medicinal Products for Human and Veterinary Use - Annex 1: Manufacture of Sterile Medicinal Products (corrected version 1. March 2010). In: Commission E, ed. Brussels2008:16.
2. European Commission. EudraLex - Volume 4 - Good Manufacturing Practice (GMP) guidelines - Chapter 1 Pharmaceutical Quality System. In: Vol SANCO/AM/sl/ddg1.d.6(2012)8603622013.
3. U.S. Food & Drug Administration (FDA). Guidance for Industry - Sterile Drug Products Produced by Aseptic Processing — Current Good Manufacturing Practice. In: FDA, CDER, CBER, ORA, eds. Rockville: FDA; 2004.
4. Task Force on Sterile Pharmaceutical Products by Aseptic Processing. Guidance on the Manufacture of Sterile Pharmaceutical Products by Aseptic Processing (1st revision of original version). In. Japan2006.
5. World Health Organization (WHO). WHO Technical Report No. 986 - Annex 2 -WHO good manufacturing practices for pharmaceutical products: main principles. In:2014:60.
6. World Health Organization (WHO). WHO Technical Report No. 996 - Annex 6 -WHO Good trade and distribution practices for pharmaceutical starting materials. In:2016:16.
7. International Council for Harmonisation (ICH). ICH Q9 - Quality risk management. In:2005.
8. International Council for Harmonisation (ICH). Pharmaceutical Development Q8(R2). In: International Conference on Harmonisation of Technical Requirements for Registration of Pharmaceuticals for Human Use, ed2009.
9. Reinmüller B, Ljungqvist B. Modern Cleanroom Clothing Systems: People as a Contamination Source. *PDA Journal of Pharmaceutical Science and Technology*. 2003;57(2):114-125.
10. Moschner C. Cleanroom Undergarments. *Cleanroom Technology*. 2002.
11. Pharmaceutical Inspection Convention, Pharmaceutical Inspection Co-Operation Scheme. *Guide To Good Manufacturing Practice For Medicinal Products - Annexes*. In. *Guide To Good Manufacturing Practice For Medicinal Products*. Geneva: PIC/S Secretariat; 2017.
12. European Medicine Agency (EMA). Annex 1 - Manufacture of Sterile Medicinal Products. In: *European Medicine Agency (EMA)*, ed2017.
13. International Organization for Standardization. ISO 14644-5: 2004 - Cleanrooms and associated controlled environments - Part 5: Operations. In. Geneva2004.
14. Institute of Environmental Sciences and Technology (IEST). IEST-RP-CC003.4:2011 - Garment system considerations in cleanrooms and other controlled environments. In. Rolling Meadows, Illinois: IEST; 2011.
15. International Organization for Standardization. ISO 13408-1:2008 - Aseptic processing of health care products — Part 1: General requirements. In. Geneva2008.
16. Institute of Environmental Sciences and Technology (IEST). IEST-RP-CC003.4:2011 - Garment system considerations in cleanrooms and other controlled environments - B2.4 Particle Dispersion Test (Body Box) Test Method. In. Rolling Meadows, Illinois: IEST; 2011.

17. Hu SC, Shiue A. Validation and application of the personnel factor for the garment used in cleanrooms. *Data Brief*. 2016;6:750-757.
18. Institute of Environmental Sciences and Technology (IEST). IEST-RP-CC003.4:2011 - Garment system considerations in cleanrooms and other controlled environments - B2.5 Helmke Drum Test Method. In: *Rolling Meadows, Illinois: IEST; 2011*.
19. European Commission. EudraLex - Volume 4 - EU Guidelines for Good Manufacturing Practice for Medicinal Products for Human and Veterinary Use - Annex 15: Qualification and Validation. In: *DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY, ed. Brussels 2015*.
20. International Organization for Standardization. ISO 11607-1:2019 Packaging for terminally sterilized medical devices -- Part 1: Requirements for materials, sterile barrier systems and packaging systems. In: *Geneva: ISO; 2019*.
21. International Electrotechnical Commission (IEC). IEC 62366-1:2015 Medical devices -- Part 1: Application of usability engineering to medical devices. In: *International Organization for Standardization*.
22. United States Pharmacopoeia (USP/NF). USP 41 <1116> Microbiological Control and Monitoring of Aseptic Processing Environments. In: *USP, ed.*

Originally published in the Journal of GXP Compliance, Volume 23, Issue 4, Jul 2019 which can be viewed at: <http://www.ivtnetwork.com>

Milenko Pavicic, PhD, started his career as a microbiologist in the pharmaceutical industry in 1994. He occupied several functions in R&D, QC and production. His specializations are in the field of contamination control related to aseptic and sterile production, microbiological quality control and (practical) training. In 2003 he founded Pavicic Pharmaceutical Microbiology (PPM). PPM supports pharmaceutical companies and hospital pharmacies in an advisory role, or by participating in or managing projects. PPM is also specialized in practical training, courses and development of tailor made e-learning modules in the field of aseptic and sterile pharmaceutical production.

Thierry Wagner has spent over 30 years working for DuPont in its polyester films and nonwovens businesses and from the 1st of October 2019 as Global Director, Regulatory & Standards – Healthcare, DuPont – Safety & Construction. He is convener of ISO TC198/WG7 “Sterilization of Health Care Products—Packaging”, chairman of the Sterile Barrier Association (SBA), member of the Parenteral Drug Association (PDA) and actively involved in various ISO and CEN technical committees on medical and pharmaceutical packaging like CEN TC102 “Sterilizers for Medical Purposes—Packaging” and ISO TC76 “Transfusion, infusion and injection equipment for medical and pharmaceutical use”. Thierry is also a member of ISO/TC 210 in charge of ISO 13485 and medical device symbols, ASTM Committee F02 and of the CEN Advisory Board for Healthcare Standards-Europe (CEN ABHS). Thierry Wagner earned a master’s degree in mechanical and process engineering from ETH Zürich in Switzerland. He is a featured speaker at international conferences and seminars on medical and pharmaceutical packaging regulatory aspects.

Annex 1 revision— Contamination risk factors

Advantages of Tyvek® IsoClean® single-use garments over reusable garments

Manufacturing of the garments

For DuPont™ Tyvek® IsoClean®, the production of raw materials is well controlled by DuPont as the sole manufacturer. For reusable garments, the risk of garment failure varies considerably across multiple manufacturers within that supply chain.

Packaging & transport of the garments

Tyvek® IsoClean® garments are packaged in a pouch with a five-year shelf life, then in double poly case liners, and finally in a single-use sealed box, significantly reducing the risk of contamination. Most reusable garments are packed in a bag with a six-month shelf life and then in a reusable crate.

Working in the cleanroom

Reusable polyester garments have a higher risk of contamination in the cleanroom due to their high permeability. Additionally, the frequent washing and sterilization of reusable garments leads to damaged fabric and larger pores, thereby increasing the release of particles and fibers into the environment. Tyvek® IsoClean® provides a consistent barrier against dry particles and microorganisms, delivering the ideal balance of protection, durability and comfort.

Washing & sterilizing the garments

Tyvek® IsoClean® garments are only washed and sterilized prior to wearing, so there is no risk of ineffective washing or sterilizing of a garment that has already been worn. Reusable garments are washed and sterilized many times over their lifetime, increasing the risk of contamination.

Repairing contaminated garments

Tyvek® IsoClean® garments are only worn once, so there is no risk of inadequately repairing damaged garments. Inadequate repairs on reusable garments can lead to holes and increased permeability. Requalifying repaired garments can also be expensive and disruptive.

Sterile single-use garment & accessory options*

Coverall

○ IC253BWHXX0025CS

Bound seams
Bound neck
Dolman sleeve
Elastic wrists
Elastic ankles
Zipper closure
25/cs
SM-4X

Bouffant

○ IC729SWH000250CS

Serged seams
Elastic headband
21.5" diameter
250/cs
Universal sizing (00)

Sleeves

○ IC501BWH000100CS

Bound seams
Covered elastic wrists, bicep
18" length
100/cs
Universal sizing (00)

Coverall

○ IC105SWHXX0025CS


Serged seams
Standard hood
Elastic hood opening
Set sleeve
Elastic wrists
Elastic ankles
Attached thumb loops
Zipper closure
Attached boots with PVC soles
25/cs
MD-4X

Hood/mask

◐ IC669BWH0001000S

Integrated hood/mask combination
Bound seams
Bound head opening
Ties with loops for fit
White hood
Blue face mask
Pleated polyethylene outer
7" wide mask
Individually packaged
100/cs
Universal sizing (00)

Boot cover

○ IC447SWHXX0100CS

Serged seams
Elastic opening
Elastic ankles
Gripper™ sole
18" high
100/cs
MD-2X

Coverall

● CC252BBUXX00250S


Bound seams
Bound neck with loop at center back
Dolman sleeve
Elastic wrists
Elastic ankles
Zipper closure
25/cs
SM-5X

Hood

○ IC668BWH000100CS

Bound seams
Full face opening
Bound hood opening
Ties with loops for fit
100/cs
Universal sizing (00)

Boot cover

○ IC457SWHXX01000S

Serged seams
Covered elastic opening
Ties at ankles
PVC sole
18" high
100/cs
SM-XL

Frock

○ IC264SWHXX0030CS

Serged seams
Bound neck
Raglan sleeve
Elastic wrists
Zipper closure
A-line
30/cs
SM-4X

Mask

○ ML7360WH0002500S

9" size
Bound Tyvek® ties
Pleated
Rayon outer facing
Metal nose piece
250/cs
Universal sizing (00)

Boot cover

○ IC458BWHXX0100CS

Bound seams
Covered elastic opening
Ties at ankles
Gripper™ sole
18" high
100/cs
SM-XL

*More styles available.
Visit safespec.dupont.com for the full portfolio.

DuPont Personal Protection
Customer service 1 800 931 3456
safespec.dupont.com
personalprotection.dupont.com
controlledenvironments.dupont.com

A global study of the performance of cleanroom garments over their life cycle

Matheus Barbosa, Jean-François Teneul

Abstract

Sterile garments for cleanroom use often present a variable performance over their entire life cycle as they are vulnerable to damage from laundering and sterilization methods. A study was conducted to understand how reusable garments perform when subjected to multiple laundering and irradiation cycles – tear strength, particle shedding, permeability, etc.. The study enables a cost comparison with single use garments.

Introduction and key concepts

In the context of a global business with ever-increasing quality standards and effectiveness requirements as well as the latest draft review of the Good Manufacturing Practice (GMP) Annex 1 of December 2017, one of the main technical areas explored is related to studying and mitigating potential process/product contamination risks against biological bodies, particulates and pyrogen agents. One key strategy for cleanroom supervisors should be to carry out complete Risk Assessments in order to map, classify, and then reduce contamination risks.

Humans are the main source of potential contamination inside cleanrooms (more than 70%), as shown by several past global studies (Akers, J. *et al.*, 2004; Ramstorp M, 2000 and Whyte and Hejab, 2007). Hence, cleanroom garments serve as the last protection barrier against controlled environment contamination by the thousands of human particles (potentially carrying microorganisms) that are shed every minute. In terms of contamination risk management, it is critical to evaluate the variables related to human contamination and cleanroom garment barrier performance besides HEPA filtering, process air flow velocity and other factors. This study aims to explore several important technical aspects of cleanroom garments that should be considered when evaluating contamination risks.

The process of wearing, laundering and sterilizing reusable cleanroom garments can impact their physical

properties and change their functionality. Laundering and wearing abrades garment fibers. Simultaneously, changes to the polymers that make up the garments can occur at the molecular level. Although routine visual inspection is often part of garment quality evaluation programmes, non-visible properties also change with time.

When selecting reusable garments for use in cleanroom environments, it is important to understand how they will perform over their intended life cycle. Consideration of all the degradation aspects should be part of the decision process for when to take reusable garments out of service, or alternatively to change to a single-use garment system. Several factors should also be considered when evaluating intrinsic risks generated by cleanroom garments, such as: particle shedding, biological/particle barrier, worker comfort and protection, durability, packaging, sterilization continuous validation – besides process and supply factors: logistics chain reliability, damages and repairs, shrinking and ergonomic fit, among others.

Physical property data are often available for new cleanroom garments; however, there are less physical property data for the remainder of the garment life cycle. To aid in garment choice, DuPont conducted a study, led by Jennifer Galvin PhD, DuPont Principal Investigator, of the physical properties of reusable cleanroom garments after a set number of laundering and gamma radiation exposure (sterilization) cycles.

Methodology

Two sets of commercially branded, reusable coveralls were purchased for testing and designated as Garment A and Garment B. Garments were made of woven polyester with integral carbon fiber for electrostatic decay properties. Garments were laundered under standard industrial settings and subsequently exposed to gamma radiation; this was considered one cycle. Garments were removed for testing after

pre-determined numbers of cycles until a total of 30 cycles had been completed.

Not all properties were tested at the same frequency. Initial properties of the garments were either measured on “as-received” garments or garments that had been laundered one time, but not exposed to gamma radiation. Parameters for garment laundering and gamma exposure were consistent throughout the study.

Garments were not worn or exposed to simulated work scenarios between cycles and the effect of routine garment “wear and tear” was not part of this study.

A summary of the garment testing methods is shown in Table I, according to IEST (Institute of Environmental Sciences and Technology), ASTM (American Society for Testing and Materials) and AATCC (American Association of Textile Colorists and Chemists) standards. Most of testing was done at third-party laboratories. Results for property testing are shown with the average and the Bonferroni confidence interval on the mean. Changes in both absolute performance and variability within the garment population may factor into formulation of end-of-life criteria.

Table 1: Test Method Summary

| Test | Test Method |
|-----------------------------------|-----------------|
| Particle shedding via Helmke Drum | IEST RP-CC003.4 |
| Particle dispersion (Body Box) | IEST RP-CC003.4 |
| Hydrostatic head | AATCC TM127 |
| Trapezoidal tear strength | ASTM D5587 |

Results and discussion

3.1 Studied parameters

Based on all experimental parameters listed above, a number of results were obtained. The object was to analyze critical limits of cleanroom garment performance in order to help end users

- evaluate which garment system solution to choose (single-use or reusable),
- or
- after how many cycles reusable garments become non-performing and have to be replaced.

The properties studied reflected key features including process protection, people protection (when needed), comfort and durability.

3.2 Polyester reaction mechanism after gamma exposure

The impact of gamma radiation exposure on a variety of polymers is well studied (Skiens, W. E, 1980). Although multiple reaction mechanisms can occur simultaneously, there is typically a predominating reaction type. The extent and type of each reaction depend on many factors and combinations of factors, including:

- Polymer composition (different polymers behave differently)
- Presence or absence of air during irradiation
- Crystallinity of the polymer and changes in crystallinity
- Physical configuration (e.g., fibre, film or tubing)
- Additional processing (e.g., laundering, calendaring or surface treatment)
- Presence of antioxidants or other additives in the polymers
- Cumulative radiation dose

The two primary reaction mechanisms that occur in polyester (PET) after exposure to gamma radiation are chain scission and cross linking (Potnis, S. P., Shetty, S. M., Rao, K. N, Prakash, J, 1969; Nair, P. D., Sreenivasan, K., and Jayabalan, M, 1988). Changes in the polymer makeup can result in changes to a garment's physical properties. To better understand which mechanism predominated under the conditions of this study, PET molecular weight was measured by size exclusion chromatography (SEC) using hexafluoro-isopropanol (HFIP) as the solvent.

Results for Garments A and B overlapped, so the data was grouped (Figure 1). Because the molecular weight of the PET decreased with laundering and exposure to gamma radiation, chain scission was the predominant mechanism. As garments were both laundered and exposed to gamma radiation (but not used by operators), this data includes the combined simultaneous impact of both factors.

This initial result indicates that reusable textile garments suffer degradation throughout their life cycle. As mentioned before, the polymer molecular weight decrease effect is a result of laundering and irradiation processes, but in actual use there are other additional effects. Repeatedly wearing and submitting garments to physical stress (standard operation moves), as well as transportation/manipulation, donning/doffing or the exposure of chemical/biological compounds could also impact and intensify degradation effects.


Figure 1: Number average polymer molecular weight (Daltons) for garments A and B

In order to analyze garment performance that result from the PET degradation process, other important physical properties were tested. These are categorized under: Process Protection, People Protection and Durability.

3.3 Process protection

The primary function of cleanroom garments is to protect a product or a process against contamination from humans (particle shedding and biological exposure) or from the garment itself (particle or linen shedding). To represent process protection, particle shedding was measured via the Helmke Drum method (Figure 2), and particle dispersion via the Body Box method (Figure 3).

3.3.1 Helmke Drum test

The best known measurable parameter for cleanroom garment cleanliness is particle shedding. This is tested in accordance with the Helmke Drum testing standard. Garments or fabric swatches are tumbled for determined cycles inside a rotating drum equipped with a standardized particle counter. The final measurement defines the shedding rate (particles/minute). Seeking to normalize test results, fabric swatches were tested, and swatch data can be evaluated for performance trends. The results show that particle shedding increased after 25 cycles or exposure to a cumulative mid-dose of 754kGy units of ionising radiation, but was fairly consistent until that point (Figure 4).


Figure 2: Helmke drum


Figure 3: Body box

Main feature

The key preliminary conclusion for these results indicates that cleanrooms that are sensitive to particle shedding should establish a monitoring program to better understand when their garments are no longer performing as required for cleanroom compliance. As observed above, after a determined number of laundering and sterilization cycles, reusable polyester materials will not only increase their magnitude of particle shedding, but also the variability of the same property. This effect can also generate an extra layer of unpredictability when determining contamination control standards and procedures – especially because the garment’s “performance breaking point” may vary based on specific application, physical stress and size fit for the operators.

A final point to consider is that the Helmke Drum test is a well-known and effective method to evaluate garment cleanliness and particle shedding from its material (polyester), but it does not indicate garment particle/bacterial barrier

performance against human shedding – which is the main source of contamination for cleanrooms. In that matter, complementary tests were conducted according to the Body Box method.

3.3.2 Body Box test

Body Box testing measures not only particle generation from the garment, but can also indicate its function as a particle barrier. The method is described in the same standard as the Helmke Drum (IEST-RP-CC003.4). In this test, a fully garbed trial subject conducts a series of movements inside a box supplied with HEPA-filtered air. Air in the box is sampled by a particle counter and shedding rate is reported as a function of activity as well as a total rate for all activities conducted during the test. This data also showed a shift in performance and variability after increased cycles of laundering and gamma radiation exposure (Figure 5).

Both the Helmke and Body Box data show an increase in both amount and

variability of shedding. Cleanroom operators who are particularly prone to particle shedding should consider establishing a monitoring programme to determine when garment performance no longer meets requirements. Particle sizes typically monitored in a cleanroom are too small to be visible to the naked eye, so visual inspection alone cannot indicate an increase in garment shedding. Other potential contamination factors related to the garments should also be studied and considered, such as the intrinsic abrasion effects of wear and tear and sterile packaging.

3.4 People protection

It is not uncommon to identify chemical and biological hazards in controlled environments or cleanrooms. In these instances, the garments not only need to perform as a process contamination barrier, but also serve as a PPE (Personal Protective Equipment) to guarantee the health and safety of the operators. Several applications and common activities may present a potential risk to workers in cleanrooms, among them:

- Oncology drugs compounding and manipulation (cytotoxic handling)
- HPAPI (High Potent Active Pharmaceutical Ingredient) manufacturing
- Hormones handling and production
- Activities with different Biosafety levels (bacteria or virus manipulation)
- Chemical products manipulation (solutions preparation, cleanroom sanitizing)
- Infectious residues (animal, human)

Several standards exist to certify chemical and biological protective garments, such as the ISO 16602 (for chemical risks) and ISO 16603/EN 14126 (for biological hazards).


Figure 4: Helmke Drum particle (greater than 0.5 microns) shedding of swatches


Figure 5: Body Box valuation by the sum of shedding for all activities

Figure 6: Hydrostatic head test

Properties such as permeation and repellency should also be considered for significant chemical/biological risks.

The hydrostatic head test (Figure 6) was used to evaluate fabric performance against an aqueous challenge. Fabric was subjected to a water column of increasing pressure until three drops penetrated the fabric. The data show a significant drop in performance as a function of exposure to laundering and gamma radiation (Figure 7). If garments are considered for incidental, light aqueous splash protection, understanding the degradation per cycle is important.

Many potential hazards are found in liquid form, so partial impermeability of the garment fabric should be considered and/or studied when selecting cleanroom equipment for applications with an exposure risk level. Part of the health and safety procedure could include establishing the garment's worker protection performance over its life cycle. As mentioned before, other international standards and certifications may assist companies to verify the chemical/biological barrier effectiveness for cleanroom garments, such as the CE Category III PPE certification and specific permeation data.

3.5 Durability

The length of a garment life cycle is also affected by the last performance parameter studied: durability. Garments should withstand normal wear and tear. Without adequate durability, garment breach is possible. Besides significant process contamination, extra costs could be generated as the repair of polyester garments is often a complicated and costly activity, sometimes not included within a company's budget or laundering contract service.

To understand the impact of laundering and exposure to gamma radiation on garment durability, trapezoidal tear strength was measured (Figure 8). Cross direction (CD) tear

strength is shown in Figure 9 while machine direction (MD) tear strength is shown in Figure 10. Often in woven garments, there are different constructions in the two directions, so differences in tear values between MD and CD are expected.

Testing showed that garment durability decreases with increasing cycles of laundering and exposure to

gamma radiation. Reducing potential impact from garment tearing is important, especially in cleanrooms and controlled environments where workers may have physical activities such as climbing stairs or bending to monitor or adjust equipment. In terms of contamination control procedures, the decay of mechanical resistance for cleanroom garments adds an extra layer of


Figure 7: Interval plot of Hydrohead (inches)


Figure 9: CD Trapezoidal tear strength (lbf - Pound-force)


Figure 10: MD Trapezoidal tear strength (lbf)

Figure 8: Tear strength test

complexity on establishing a monitoring standard to determine each garment's life cycle end point; variables such as workers' activities and potential external abrasion factors also need to be considered.

Conclusions

People are the main source of contamination in controlled environments, and cleanroom garments are the main and last barrier to protect critical processes and products. The data outlined here demonstrate that garment properties do change after several laundering and gamma exposure cycles. These changes are not always visible to the naked eye, so visual garment inspection alone may not be sufficient to understand garment performance. Based on these findings, the following guidelines are recommended:

- Even though tests were conducted with sterilization via gamma irradiation, several studies show that other sterilization methods also present abrasion and degradation effects over the garment's life cycle. Autoclaving, for example, uses a physical process that may degrade polyester composition after several utilization cycles (Nair, P. D., Sreenivasan, K, 1984). It is important to consider that continuous laundering (shrinking & expanding), and wearing of cleanroom garments also play a considerable part in the structure degradation of fabrics.
- Consider performance data over the entire garment life cycle. If not available, question your cleanroom garment provider or assess the risk of not having control of your garments' system in place;
- When cleanroom garments also need to perform as Personal Protective Equipment, companies should consider looking for specific technical data and certifications that would enhance worker safety and protection. Asking garment providers for permeation data for specific risks or barrier technical claims might be an effective strategy;
- Enact testing protocols to monitor the performance of garments as they age, based on the risk assessments and needs of each individual cleanroom. Parameters should not only consider particle shedding and cleanliness of the garment itself,

but also its barrier effectiveness against human contamination (particle shedding and biological filtering) and sterility validation assurance. Then, establish criteria for taking garments out of service when they no longer meet functionality requirements;


- It is also important to continuously map, evaluate and control the risk of the entire garments system value-chain: from the fabric weaving and sourcing, garment assembly, packaging and sterilization and, if applicable, the laundry process as well.

It ought to be noted that since garment requirements vary by cleanroom operation, establishing initial and ongoing fitness for use is the responsibility of the end user. Garment assessment may require evaluation of additional information beyond what is presented here. For example, seams and closures may have lower barrier properties than fabric. Properties of garments and fabrics subjected to other conditions, including different sterilization methods, may also vary.

In conclusion, when assessing risks related to potential contamination in controlled environments, cleanroom companies should question and require their suppliers to support their quality claims with continuous technical data, risk mitigation and process control.

References

- Akers, J. *et al.*, Development of a Highly Automated PAT- Compatible Aseptic Processing System, paper presented at the A3P congress, Biarritz, France, 2004.
- Galvin, J (DuPont de Nemours). A Life Cycle Assessment of Reusable Garment Properties, Cleanroom Technology, June 2017
- Skien, W. E. Radiat. Sterilizing radiation effects on selected polymers. *Phys. Chem.*, 1980, 15, p. 47-57
- Potnis, S. P., Shetty, S. M., Rao, K. N. and Prakash, J. *Die Angew.* Effect of γ -radiation on synthetic fibres – II irradiation in the presence of chemicals. *Markromol. Chemie*, 1969, 6, p. 127-135
- Nair, P. D., Morphological changes of poly(ethylene terephthalate) on multiple steam sterilization. *Clinical Materials*, 1990, 5, p. 43
- Nair, P. D., Sreenivasan, K. Effect of steam sterilization on polyethylene terephthalate. *Biomaterials*, 1984, 5, p. 305
- Nair, P. D., Sreenivasan, K., and Jayabalan, M. Multiple gamma radiation sterilization of polyester fibres. *Biomaterials*, 1988, 9, p. 335-338
- Whyte, W. and Hejab, M., Particle and microbial airborne dispersion from people, *European Journal of Parenteral and Pharmaceutical Sciences*, 12 (2). pp. 39-46. ISSN 0964-4679, 2007

Matheus Barbosa, graduated as a Chemical Engineer at the Polytechnic School of the University of São Paulo, Brazil, and has over five years of experience in DuPont as a personal protection product specialist and technical expert for controlled environments and cleanroom garments. He is currently Latin America Marketing Manager for DuPont Personal Protection covering chemical, thermal and mechanical safety equipment.
Contact: Matheus.Rodrigues.Barbosa@dupont.com, www.epi.dupont.com.br

Jean-François Teneul is Global Life Sciences Market Leader for Personal Protection at DuPont de Nemours, Luxembourg. His current role is to bring innovative cleanroom clothing solutions to ensure a safe manufacturing environment and to reduce the risk of contamination in Controlled Environments. He is experienced in Infection and Contamination Control and in Safety and Protection. He started in Health Care and Medical Devices industry and moved to Personal Protective Equipment (PPE). He served as chairman of EDANA (European Disposables and Nonwovens Association) for the MEDECO group and is Six Sigma Black Belt certified.
Contact: Jean-Francois.Teneul@dupont.com, www.ipp.dupont.com

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
DuPont™ SafeSPEC™
www.safespec.dupont.de

Tyvek IsoClean vs Reusable Polyester Textile - Droplet Tests.mp4

This slide was originally a mp4 file and could not be converted into PDF.
You can find the original file in the .ZIP that was downloaded.

If you're viewing this in Adobe Reader, this file can be opened automatically by clicking [**here**](#).


Tyvek IsoClean vs Reusable Polyester Textile - Hydrohead Tests.mp4

This slide was originally a mp4 file and could not be converted into PDF.
You can find the original file in the .ZIP that was downloaded.

If you're viewing this in Adobe Reader, this file can be opened automatically by clicking **[here](#)**.

GRÖSSENTABELLE

| Größe | Brustumfang (cm) | Körpergröße (cm) | Brustumfang (Zoll) | Körpergröße (Fuß /Zoll) |
|-------|------------------|------------------|--------------------|-------------------------|
| XXS | 68 - 76 | 150 - 158 | 27 - 30 | 4'11" - 5'2" |
| XS | 76 - 84 | 156 - 164 | 30 - 33 | 5'1" - 5'5" |
| S | 84 - 92 | 162 - 170 | 33 - 36 | 5'4" - 5'7" |
| M | 92 - 100 | 168 - 176 | 36 - 39 | 5'6" - 5'9" |
| L | 100 - 108 | 174 - 182 | 39 - 43 | 5'8" - 6'0" |
| XL | 108 - 116 | 180 - 188 | 43 - 46 | 5'11" - 6'2" |
| 2XL | 116 - 124 | 186 - 194 | 46 - 49 | 6'1" - 6'4" |
| 3XL | 124 - 132 | 192 - 200 | 49 - 52 | 6'3" - 6'7" |
| 4XL | 132 - 140 | 200 - 208 | 52 - 55 | 6'7" - 6'10" |
| 5XL | 140 - 148 | 208 - 216 | 55 - 58 | 6'10" - 7'1" |
| 6XL | 148 - 156 | 208 - 216 | 58 - 61 | 6'10" - 7'1" |
| 7XL | 156 - 162 | 208 - 216 | 61 - 64 | 6'10" - 7'1" |