

# EU Type Examination Certificate

This is to certify that:

Berner International GmbH  
Werner-von-Siemens-Str. 19  
Elmshorn  
Schleswig-Holstein  
25337  
Germany

Holds Certificate Number:


CE 715808

In respect of:

**Model cleo® saphir and cleo® saphir for Codan Chemoprotect® Protective Clothing range of gowns, oversleeves and overshoes  
To EN 14126:2003 and EN 14605:2005+A1:2009**

on the basis that BSI carried out the relevant Type Examination procedures under the requirements with the Regulation (EU) 2016/425 of the European Parliament and Council relating to Personal Protective Equipment Regulation (PPE) Annex V (Module B) and meets the relevant health and safety requirements specified in Annex II

For and on behalf of BSI, a Notified Body for the above Regulation (Notified Body Number 2797):



Drs. Dave Hagenaaars, Managing Director

First Issued: 2019-07-29

Latest Issue: 2021-02-05

Effective Date: 2021-02-05

Expiry Date: 2024-07-29

Page: 1 of 4



...making excellence a habit.™

# EU Type Examination Certificate

No. CE 715808

## Product Specification

**Range:** cleo® saphir and cleo® saphir for Codan Chemoprotect® Protective Clothing

**Product Description:** The cleo® saphir and cleo® saphir for Codan Chemoprotect® protective equipment range is manufactured from spun polypropylene fleece with a coating of polyethylene. The range comes in two types of colour and design. The protective equipment is impervious to liquids on the sleeves and front parts and it is a latex-free product. It is intended for single use only.

## Models in Range:

Model Number	Description
6700/6701	cleo® saphir gowns with knitted cuffs, blue, unsterile/sterile
6800/6801	cleo® saphir gowns with knitted cuffs, blue, unsterile/sterile
6900/6901	cleo® saphir gowns with knitted cuffs, blue, unsterile/sterile
100072/100073	cleo® saphir gowns with knitted cuffs, blue, unsterile/sterile
6000/6001	cleo® saphir oversleeves with knitted cuffs, blue, unsterile/sterile
125501/ 125502/ 125505	cleo® saphir overshoes
6500/6600	cleo® saphir gowns with elasticated cuffs, light blue, unsterile/sterile
6550/6650	cleo® saphir gowns with elasticated cuffs, light blue, unsterile/sterile
6200/6300	cleo® saphir oversleeves with elasticated cuffs, light blue, unsterile/sterile
100253/79.5200, 100254/79.5201, 100255/79.5202, 100256/79.5203, 100257/79.5204, 100258/79.5205	cleo® saphir for Codan Chemoprotect® gowns with knitted cuffs, blue, unsterile/sterile
100259/79.5500, 100260/79.5501	cleo® saphir for Codan Chemoprotect® oversleeves with knitted cuffs, blue, unsterile/sterile
100261/26.2514	cleo® saphir for Codan Chemoprotect® overshoes, blue

First Issued: 2019-07-29

Latest Issue: 2021-02-05

Effective Date: 2021-02-05

Expiry Date: 2024-07-29

Page: 2 of 4

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.

# EU Type Examination Certificate

No. CE 715808

## Product Specification (continued)

**Sizes:** Gowns - S, M, L, XL  
Oversleeves / overshoes – single size

**Classification:** Category III PPE (Complex)  
Type PB[4]-B

**Technical Specification:** EN 14126:2003 Protective clothing — Performance requirements and tests methods for protective clothing against infective agents  
EN 14605:2005+A1:2009 Protective clothing against liquid chemicals - performance requirements for clothing with liquid-tight (Type 3) or spray-tight (Type 4) connections, including items providing protection to parts of the body only (Types PB [3] and PB [4])

First Issued: 2019-07-29

Latest Issue: 2021-02-05

Effective Date: 2021-02-05

Expiry Date: 2024-07-29

Page: 3 of 4

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.

# EU Type Examination Certificate

No. CE 715808

## Certificate Administration Details

**Technical File Reference:** 'TECHNICAL DOCUMENTATION FOR PPE, cleo® saphir Protective Clothing'

## Certificate Administration Record and BSI internal Technical File Review reference

Issue date	Comments	BSI Project Ref.
February 2021	Transition from 0086 NB to 2797 NB	2797:21:3319175
July 2019	First issue	0086:19:3057686

Note: The Certificate holder is responsible for keeping the Notified Body advised of changes to any aspect of the overall process used in the manufacture of the product

## Monitoring of manufactured PPE

The EC Quality of Production by Monitoring requirements, Module C2, for the product are referenced in BSI issued Certificate Number CE 715809

First Issued: 2019-07-29

Latest Issue: 2021-02-05

Effective Date: 2021-02-05

Expiry Date: 2024-07-29

Page: 4 of 4

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.