

PRODUCT INFORMATION

Berner Laboratory Overshoes



Area of application and features

- + **Area of application:** Single use protective overshoes for the use in the laboratory or similar areas. Elasticated seam for a good fit, suitable for all shoes sizes.
- + **Overshoe replacement interval:** Single use!
- + **Before use:** Check for damage! Do not use damaged overshoes!!

Type

Size	universal
Item no.	100086
Length	~ 41 cm
Colour	blue
Elasticated seam	elasticated seam
Packaging unit	100 pairs

Material properties

Material	polyethylene
Material thickness	~ 0.02 mm
Weight (pair)	~ 6 g

Quality management system

Our **quality management system** is **tested and certified** in accordance with DIN EN ISO 9001:2015 by the TÜV Product Management Service GmbH (a certification body accredited by the German Accreditation Council). Regular **audits and production site inspections** guarantee the quality of our products.

Storage and transport conditions

- + Dark (protect from direct UV light and sunlight)
- + Cool (+5 to +40°C)
- + Dry
- + No contact with pointed and/or sharp objects

Shelf Life

- + 3 years from the date of manufacture