

PRODUCT INFORMATION

Light Inhibiting Bag

Area of Application and Features

- + **Area of application:** UV light inhibiting bags with a resealable micro-snap seal for the storage and transportation of substances which are sensitive to light, e.g. cytostatic drugs.
- + **Features:** The material is opaque and nontransparent. Resealable.

Material Properties and Versions

Dimension	150 x 200 mm (± 5 mm)	230 x 300 mm (± 5 mm)
Item no. / Packaging Unit	12 01 50 / 300 pieces	12 02 30 / 300 pieces
Item no./ Packaging Unit	12 15 02 00 / 1000 pieces	12 23 03 00 / 1000 pieces
Weight	~ 3 g	~ 9 g
Material	Polyethylene (LDPE)	
Material thickness	~ 50 µm	
Colour	opaque brown	
UV Stabilisator / Absorber	Protection against UV radiation from 280 – 400 nm (UVA & UVB), equal to light energy of at least 90 kly	

Quality Management System

Our quality management system is tested and certified by TÜV Product Service GmbH (a certification body accredited by the German Accreditation Council in accordance with DIN EN ISO 9001:2015. Regular audits and production site inspections guarantee the quality of our products.

Storage and Transport Conditions

- + Cool (not above 25°C)
- + Dry (relative humidity 30% - 60%)
- + No contact with pointed and/or sharp objects

Shelf Life

2 years from the date of manufacture

Distributor

Berner International GmbH, Werner-von-Siemens-Str. 19, 25337 Elmshorn, Germany

Tel: +49 4121 43560, Fax: +49 4121 435620

info@berner-safety.de, www.berner-safety.de