

# PRODUCT INFORMATION

## Aluminum cleaning tool telescopic

### Area of application and features

- + Device for efficient and straight-forward cleaning of isolators and safety cabinets, perfect for hard to reach places. For simple cleaning of wall, floor and ceilings.
- + Ergonomic, maneuverable aluminum tool.
- + Handle grip for ergonomic cleaning.
- + The cleaning tool can be sterilized autoclavable 121 °C.
- + Suitable for sterile cleanrooms.
- + Quick-change clip system.
- + It contains two parts: telescopic mop handle and mop frame.

### Material

<b>Material</b>	Aluminum and Polyoxidemethylene
<b>Color</b>	Gray and white
<b>Length</b>	min. 55,2 cm – max. 88,5 cm
<b>Width</b>	19,0 cm
<b>Diameter</b>	4,5 cm
<b>Weight</b>	540 g
<b>Packaging</b>	1 Set / Bag



### Product information

<b>Item No.</b>	204882
<b>Sterilization</b>	Autoclavable at 121 °C for 20 min or 134 °C for 5 min.
<b>Cleanroom classification</b>	ISO 4-8 GMP A/B, C, D
<b>Resistance to chemicals</b>	Suitable for use with common biocides, detergents and alcohol. Not suitable for use with sodium hydroxide and hydrochloric acid.