

PRODUCT INFORMATION

Stainless-steel cleaning tool telescopic

Area of application and features

- + Device for efficient and straight-forward cleaning of isolators and safety cabinets, perfect for hard to reach places. For simple cleaning of wall, floor and ceilings.
- + Ergonomic, manoeuvrable stainless-steel tool, durable.
- + Handle grip for ergonomic cleaning.
- + The cleaning tool can be sterilized autoclavable 121 °C.
- + Suitable for sterile cleanrooms.
- + Quick-change clip system.
- + It contains two parts: telescopic mop handle and mop frame.

Material

Material	Stainless-steel
Color	Gray and white
Length	min. 55,0 cm – max. 89,0 cm
Width	19,0 cm
Diameter	4,5 cm
Weight	670 g
Packaging	1 Set / Bag



Product information

Item No.	204883
Sterilization	Autoclavable at 121 °C for 20 min or 134 °C for 5 min.
Cleanroom classification	ISO 4-8 GMP A/B, C, D
Resistance to chemicals	Suitable for use with common biocides, detergents and alcohol.