

PRODUCT INFORMATION

Cleaning tool for Isolators & Safety Cabinets

Applications and Features

- + Device for efficient and straight-forward cleaning of isolators and safety cabinets, perfect for hard to reach places. For simple cleaning of wall, floor and ceilings.
- + Ideally suited for use with BERNER Isysoft cleaning wipes (available in sterile, non-sterile or latex-free).
- + Quick and easy fastening or replacing of wipes using 4 polyethylene clips.
- + With silicone pad for effective cleaning even on uneven surfaces.
- + The cleaning tool can be sterilised* (gamma radiation, fumigation with hydrogen peroxide / ethylene oxide, etc.)* and is therefore suitable for use in cleanroom classes A and B.
*not autoclavable

Material Properties Cleaning Tool I

Total length	ca. 460 mm
Size	ca. 170 x 90 mm
Length handle	Ca. 430 mm



To order

	Order No.	
Cleaning tool I	802645	
Isysoft cleaning cloths, available separately:	Non-sterile	Sterile*
Isysoft (4 x 25 wipes) / (50 x 2 wipes)*	121500	121501
Isysoft latexfree (4 x 25 wipes) / (50 x 2 wipes)*	121507	121508
Cleanroom 100 (150 wipes) / (25 wipes)*	121740	121471
Extra clip (1 piece)	802824	