



**CLEO® SAPHIR - PROTECTIVE CLOTHING**  
Highest protection & reliability in the laboratory



## Product highlights

### STRONG LAMINATE AND OPTIMUM SEAM TECHNOLOGY

Our proven Z+ B+ protective clothing, including gowns, oversleeves and overshoes, has been improved and will be sold under the new name of Cleo® saphir. Our Cleo® saphir protective clothing was also subject to a new type test examination, so that our partial body protection products are now certified according to the latest EU-regulation for personal protective equipment (EU) 2016/425.

- PPE of the highest category III**  
 Type PB-[4]/ PB-[4]-B to protect against chemicals, cytostatics and biological agents
- Best infection protection**  
 Maximum resistance against biological agents such as viruses, bacteria and fungi
- Strong coating**  
 Highest protection level against all tested chemicals and cytotoxic drugs
- First class manufacturing conditions**  
 Production in a controlled environment and preparation for sterilisation in ISO class 7
- Certified according to the latest EU-regulation for personal protective equipment (EU) 2016/425**
- Optimum seam technology**  
 of the sealed sleeve seam, tested with critical test substance Carmustine, guarantees the highest protection for at least 8 hours



## Protective gown

Proven blue partial body protection with knitted cuffs, a completely closed front, extra long cut and a raised collar. Improved 3-fold ultrasonic sealed seams on the sleeves for greater protection. The combination of liquid-impermeable and breathable materials ensures comfortable and safe working.

**To order**  
cleo® saphir protective gown with knitted cuffs

Size	S	M	L	XL
Order No. (non-sterile / 15 pcs)	6700	6800	6900	
Order No. (non-sterile / 10 pcs)				100072
Order No. (sterile / 10 pcs)	6701	6801	6901	100073



## Protective gown

Proven light blue partial body protection with elasticated cuffs, a completely closed front, extra long cut and a raised collar. Improved 3-fold ultrasonic sealed seams on the sleeves for greater protection. The combination of liquid-impermeable and breathable materials ensures comfortable and safe working.

**To order**  
cleo® saphir protective gown with elasticated cuffs

Size	M	L
Order No. (non-sterile / 15 pcs)	6500	6550
Order No. (sterile / 10 pcs)	6600	6650





## Protective oversleeves

Proven blue partial body protection with knitted cuffs and cone-shaped cut for wearing comfort. Improved 3-fold ultrasonic sealed seams for greater protection. The liquid-impermeable materials ensure safe working.

**To order**

cleo® saphir protective oversleeves with knitted cuffs

Size	uni
Order No. (non-sterile / 50 pairs)	6000
Order No. (sterile / 40 pairs)	6001



## Protective oversleeves

Proven light blue partial body protection with elasticated cuffs and cone-shaped cut for more wearing comfort. Improved 3-fold ultrasonic sealed seams for greater protection. The liquid-impermeable materials ensure safe working.

**To order**

cleo® saphir protective oversleeves with elasticated cuffs

Size	uni
Order No. (non-sterile / 50 pairs)	6200
Order No. (sterile / 30 pairs)	6300



## Protective overshoes

Proven blue partial body protection with slip-resistant sole. The universal fit and the elasticated ankles guarantee a comfortable fit and are easy to don and doff.

**To order**

cleo® saphir protective shoes with elasticated ankles

Size	uni
Order No. (non-sterile / 40 pairs)	125501
Order No. (sterile / 30 pairs)	125502



Berner International GmbH  
Werner-von-Siemens-Str. 19  
25337 Elmshorn  
Germany

Tel. +49 4121 4356-0  
Fax +49 4121 4356-20  
Email [info@berner-safety.de](mailto:info@berner-safety.de)

[www.berner-safety.de](http://www.berner-safety.de)